

Steps to Capturing Patient Wristband:

Step #1

The Patient Wristband button is located in the narrative tab of the EPCR.

The screenshot shows the EPCR Narrative tab interface. At the top, there is a navigation bar with icons for Incident, Patient, Disposition, Vitals, Procedures, Med Admin, Narrative, and Report. The Narrative tab is selected. Below the navigation bar, there is a large text area labeled "TEST NARRATIVE." Below this, there are several input fields and buttons. The "Patient Found" dropdown is set to "Lying". The "Patient Transport" dropdown is set to "W/O INCIDENT OR ADDITIONAL CHANGE IN CONDITION". The "Author Signature" field is empty. The "Hosp. Team Activation" dropdown is set to "No". The "Hosp." dropdown is set to "Summa Barberton". The "Hosp. Disposition" dropdown is set to "In Hosp Distination". The "Patient Wristband" button is highlighted with a callout bubble labeled "Step #3".

Step #2

Select the Hospital destination


Step #3

Tap the Patient Wristband button.



Step #4

The camera screen will open and will automatically capture barcode once it is able to get into focus.



IMS will only accept the patient barcode into the system. Scanning any other barcode will result in an error message.

Step #5

Complete documentation and sync as normal.

